

Mint Top Tube (preferred) -OR- Gold Top Tube		Gold Top Tube ONLY	Purple Top Tube ONLY	Pink Top Tube ONLY	Blue Top Tube ***Platelet Function Testing: Lab Draw Only, Call Lab
4 ml lithium heparin -OR- 5 ml gold		5 ml gold	3 ml purple	6 ml pink	2.7 ml (90% minimum fill) blue
Amylase/Lipase	Phosphorus	AFP	CBC	ABO RH	APTT
Aceta	PTN	CA 15-3	H&H	TYPE AND SCREEN	Fibrinogen
Alcohol	Renal Panel	CA 19-9	PLATELET ONLY	Type & Crossmatch	PT INR
Bascimeta	Salicylate	CA 27.29	Sed Rate (ESR)	DAT	D DIMER HS
Bili Group	Theo	CA 125	RETIC	TRXN	FDP
BUN/Creat	Tobra	CEA	Any Combination of the above tests can go in the same tube, Others need separate tubes	Draw 1 pink top and 1 red top tube	HIT Antibody
Carba	TSH, T4, FT3, FT4	Cortisol			Anti-Xa
Compmeta	Tup	FSH			DIC SCREEN -Draw 2 Blue top tubes (90% min fill) and 1 Purple top tube
Digoxin	Val Acid	Haptoglobin			
Ferritin	Vanco	Hepatitis	ABORH Recheck		Platelet ARU (Aspirin) Test*** Draw 2 Blue tops with a discard
Folate	Vitamin B12	HIV	CORD		- No Nurse Draws -
Gentamycin	Vitamin D	Iron Group(IRON, TIBC)	Fet Bleed		Platelet PRU (Plavix) Test*** Draw 2 Blue tops with a discard
GTT		Ketone	Fet Hgb		- No Nurse Draws -
HCG		LH	SICKLE SCREEN		
HCY		Lithium			
IGA		PCT	BNP (B Type Pep)		
LDH		Prealbumin	HGBA1C		
Lipid Panel		Prolactin	IPTH		ON ICE (Bring to Lab Immediately)
Liver Panel		PSA			
Magnesium		Rubella			Ammonia Purple 3 ml
Phenobarb		Serology Tests			
		Serum Osmo			
		T3	Red Top		
4 ml lithium heparin ONLY		Testosterone	RPR		Lactate 2 ml grey top
High Sensitivity Troponin		Tryptase	Keppra		

REV 1/31/24 CT

Mint Top Tube (preferred) -OR- Gold Top Tube		Gold Top Tube ONLY	Purple Top Tube ONLY	Pink Top Tube ONLY	Blue Top Tube ***Platelet Function Testing: Lab Draw Only, Call Lab
4 ml lithium heparin -OR- 5 ml gold		5 ml gold	3 ml purple	6 ml pink	2.7 ml (90% minimum fill) blue
Amylase/Lipase	Phosphorus	AFP	CBC	ABO RH	APTT
Aceta	PTN	CA 15-3	H&H	TYPE AND SCREEN	Fibrinogen
Alcohol	Renal Panel	CA 19-9	PLATELET ONLY	Type & Crossmatch	PT INR
Basicmeta	Salicylate	CA 27.29	Sed Rate (ESR)	DAT	D DIMER HS
Bili Group	Theo	CA 125	RETIC	TRXN	FDP
BUN/Creat	Tobra	CEA	Any Combination of the above tests can go in the same tube, Others need separate tubes	Draw 1 pink top and 1 red top tube	HIT Antibody
Carba	TSH, T4, FT3, FT4	Cortisol			Anti-Xa
Compmeta	Tup	FSH			DIC SCREEN -Draw 2 Blue top tubes (90% min fill) and 1 Purple top tube
Digoxin	Val Acid	Haptoglobin			
Ferritin	Vanco	Hepatitis	ABO Recheck		Platelet PRU (Plavix) Test*** Draw 2 Blue tops with a discard
Folate	Vitamin B12	HIV	CORD		- No Nurse Draws -
Gentamycin	Vitamin D	Iron Group(IRON/TIBC)	Fet Bleed		Platelet ARU (Aspirin) Test*** Draw 2 Blue tops with a discard
GTT		Ketone	Fet Hgb		- No Nurse Draws -
HCG		LH	Sickle Screen		
HCY		Lithium			
IGA		PCT	BNP (B Type Pep)		
LDH		Prealbumin	HGBA1C		
Lipid Panel		Prolactin	IPTH		ON ICE (Bring to Lab Immediately)
Liver Panel		PSA			
Magnesium		Rubella			Ammonia Purple 3 ml
Phenobarb		Serology Tests			
		Serum Osmo			
		T3	Red Top		
4 ml lithium heparin ONLY		Testosterone	RPR		Lactate 2 ml grey top
High Sensitivity Troponin		Tryptase	Keppra		

REV 1/31/24 CT